


1. If you don't see a palette, click [Workbook] in the menu



2. Click to display or remove whichever palette you like.


3. When you are ready to take and display data, select [Table & Graph]




Clicking [**Hardware Setup**] should give an image of the PASCO 850 Universal Interface.

If not, click [**Choose Interface**].

You might need to open the [**Manually Choose**], then check **PASCO 850 Universal Interface**.



The PASCO Interface will allow different measuring devices. To plug in the motion sensor, **right-click** on the left-most digital input, then scroll down and select the [**Motion Sensor II**].



Untitled.cap\*

File Edit Workbook Display Journal Help

Page #1

Tools

- Hardware Setup
- Data Summary
- Calibration
- Calculator
- Signal Generator

[Table title here]

	Time (s)	Position (m)
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		
11		
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		
23		
24		
25		
26		
27		
28		

<No Data Selected>

Position (m)

<Select Measurement>

- Add Similar Measurement
- QuickCalc
- Motion Sensor II
- Position (m)
- Velocity (m/s)
- Acceleration (m/s<sup>2</sup>)
- Equations/Constants
- Constants
- Time
- Time (s)
- Index
- Index

Click each <Select Measurement>  
In the pop-up menu, select the appropriate variable.

[Graph title here]

Time (s)

0 2 4 6 8 10

0 5

Displays

- Graph
- Scope
- FFT
- Histogram
- Digits
- Meter
- Table
- Text Box
- Text Entry Box
- Image
- Movie


Controls

Record Continuous Mode 00:00.00 Ready

Motion Sensor II 20.00 Hz

Recording Conditions

Delete Last Run




When you are ready, click [Record] to begin collecting data.


The [Record] button will become a [Stop] button while the data is being collected.

The graphs and table should fill with data as it is collected.

When you are done, click [Stop] to end the data collection.


















[Table title here]

		Run #1	Run #1
	Time (s)	Position (m)	Velocity (m/s)
81	4.000	0.660	0.21
82	4.050	0.670	0.21
83	4.100	0.681	0.22
84	4.150	0.692	0.22
85	4.200	0.703	0.22
86	4.250	0.714	0.23
87	4.300	0.725	0.23
88	4.350	0.737	0.23
89	4.400	0.748	0.23
90	4.450	0.760	0.23
91	4.500	0.771	0.24
92	4.550	0.783	0.24
93	4.600	0.795	0.24
94	4.650	0.807	0.24
95	4.700		
96	4.750		
97	4.800		
98	4.850		
99	4.900		
100	4.950		
101	5.000		
102	5.050		
103	5.100		
104	5.150		
105	5.200		
106	5.250		
107	5.300	0.974	0.27
108	5.350	0.988	0.28

When you click here to create a best-fit, you get a menu to select the type of curve to fit the data with. You will usually want either linear or quadratic fit.



- Proportional: At
- Linear:  $mt + b$
- Weighted Linear:  $mt + b$
- Quadratic:  $At^2 + Bt + C$
- Cubic:  $A + Bt + Ct^2 + Dt^3$
- Polynomial:  $A + Bt + Ct^2 + \dots + Gt^6$
- Power:  $A(t-t_0)^n + B$
- Inverse (no offset):  $A/t + B$
- Inverse:  $A/(t-t_0) + B$
- Inverse Square (no offset):  $A/t^2 + B$
- Inverse Square:  $A/(t-t_0)^2 + B$
- Inverse Power:  $A/(t-t_0)^n + B$
- Natural Exponential:  $Ae^{(-Bt)} + y_0$
- Natural Logarithm:  $A \ln(B(t-t_0)) + C$
- Base-10 Exponential:  $A 10^{(Bt)} + C$
- Base-10 Logarithm:  $A \log(B(t-t_0)) + C$
- Inverse Exponent:  $A(1 - e^{(-B(t-t_0))}) + C$
- Sine:  $A \sin(\omega t + \phi) + C$
- Sine Series:  $A_1 \sin(\omega_1 t - \phi_1) + A_2 \sin(\omega_2 t - \phi_2) + \dots + B$
- Damped Sine:  $Ae^{(-Bt)} (\sin(\omega t + \phi)) + C$
- Cosine Squared:  $A \cos^2(Bt + \phi) + C$
- Gaussian:  $A e^{-(t-t_0)^2 / (2B^2)} + y_0$
- Normalized Gaussian:  $(1/(\sigma \sqrt{2\pi})) e^{-(t-t_0)^2 / (2\sigma^2)} + y_0$
- User Defined:  $f(t)$



Velocity versus Time (Constant Acceleration)

