

Physics Lab Set-up

Thomas More College

Physics Lab Set-up

Thomas More College

Joe Christensen
Thomas More College

Credit for MathBookXML / PreTeXt format:

Robert A. Beezer

Latest update: October 12, 2017

Edition: Annual Edition 2017

Website: [TMC Physics](#)

© 2017–2018 J. Christensen

Permission is granted to copy, distribute and/or modify this document under the terms of the GNU Free Documentation License, Version 1.2 or any later version published by the Free Software Foundation; with no Invariant Sections, no Front-Cover Texts, and no Back-Cover Texts. A copy of the license is included in the appendix entitled “GNU Free Documentation License.” All trademarks™ are the registered® marks of their respective owners.

Acknowledgements

I would like to acknowledge the following reviewers for their helpful comments and suggestions.

- Physics adjunct instructors who provided feedback about the labs
 - Tom Neal, Physics Adjunct
- Physics faculty who designed and contributed to the design of the lab experiments
 - Jack Wells
 - Dr. Wes Ryle
 - Dr. Jeremy Huber

Download the PDF here

A PDF of this document can be found at <http://physics.thomasmore.edu/PHY121Lab/TMC-lab-setup.pdf>. (288 kB)

Contents

Acknowledgements	v
Download the PDF here	vii
1 PHY 121L: General Physics (algebra-based, fall)	1
1.1 Meaningful Measurements	1
1.2 Standard Deviation	2
1.3 Constant Acceleration	3
1.4 Newton's 2 nd Law on a Linear Track with the Sonic Ranger	4
1.5 Dry Sliding Friction	5
1.6 Centripetal Force	8
1.7 Hooke's Law / Springs / Peer Review	9
1.8 Pendulum / Peer Review	12
1.9 Next Lab	14
1.10 Human Forearm	14
2 PHY122L: General Physics (algebra-based, spring)	17

Class 1

PHY 121L: General Physics (algebra-based, fall)

1.1 Meaningful Measurements

Location	Equipment	Notes
For Each Lab Station		
AA13	1 Vernier caliper (Figure 1.1.2)	Ask faculty if they want digit or analog before and after lab, verify digital calipers are turned off
AA14	1 micrometer (Figure 1.1.2)	before and after lab verify jaws are not tight
AF34-14	1 metric ruler	could be a 1-foot ruler or a 0.25-meter stick
AF55-19	1 tall, skinny, graduated cylinder	at least 2 of the 3 objects should fit inside the graduated cylinder
AE21	string	sufficient string to tie onto the objects in order to immerse them in the cylinder and retrieve them
AE82	3 objects to measure	Ask faculty which objects they want (smooth sphere, rough sphere, cube, block of metal, irregular shape, etc)
At the front for students to share		
S224	at least 1 digital scale	the available scale(s) should be able to weigh the chosen objects
AF36-4		

Table 1.1.1: Equipment Needed: Meaningful Measurements




Figure 1.1.2: The Vernier Caliper




Figure 1.1.3: The Micrometer Caliper

(Updated: October 12, 2017)

A digital version of the lab should be found at <http://physics.thomasmore.edu/PHY121Lab/c-meaningful-measurements.html>

A PDF version of the write-up might be found at [Measurement.pdf \(291 kB\)](#)

1.2 Standard Deviation

Location	Equipment	Notes
At the front for students to share		
AL32-5	several boxes/bags of pennies	in sufficient number for each student to have up to 50 pennies

Table 1.2.1: Equipment Needed: Standard Deviation

(Updated: September 6, 2017)

A digital version of the lab should be found at <http://physics.thomasmore.edu/PHY121Lab/c-standard-deviation.html>

A PDF version of the write-up might be found at [StDev.pdf \(232 kB\)](#)

1.3 Constant Acceleration

Location	Equipment	Notes
For Each Lab Station		
S224	1 motion sensor (same as “sonic ranger”)	Should have a black-yellow plug
AF12 (tube)	1 track	.
AF15	Wood Squares	probably 2-3, used to prop up one end of track
AF22-2-13	1 cart with “sail”	these are in a large box labelled “DYNAMIC CARTS AF22-2-13”
AF34-14	ruler	used to level the track
either AF35 (shelf) or AE82 (drawer)	1 metal ball (any size)	used to level the track
-	Pasco	Computer
At the front for students to share		
AA41 or AA42 (drawers)	1 gravity protractor	This is the large yellow protractor
AF44	1 pendulum bob	

Table 1.3.1: Equipment Needed: [Constant Acceleration](#)

(Updated: October 12, 2017)

A digital version of the lab should be found at <http://physics.thomasmore.edu/PHY121Lab/c-acceleration.html>

A PDF version might be found at [Acceleration.pdf](#)

1.4 Newton's 2nd Law on a Linear Track with the Sonic Ranger

Location	Equipment	Notes
For Each Lab Station		
S224	1 motion sensor (AKA "sonic ranger")	Should have a black-yellow plug
AF12 (tube)	1 track	.
AF15	2-3 small wooden squares	used to level the track
AF22-2-13	1 cart	these are in a large box labelled "DYNAMIC CARTS AF22-2-13"
AF22-2-13	1 wooden cart-block	these are in a large box labelled "DYNAMIC CARTS AF22-2-13"
AF22-2-13	light plastic bucket	These might already be attached to the string in a large box labelled "DYNAMIC CARTS AF22-2-13"
(attached to plastic bucket?)	string	There should be pre-cut string that is long enough to reach from the cart, over the pulley and to a hanging mass. About one meter long. This is likely in a large box labelled "DYNAMIC CARTS AF22-2-13"
AF22-2-13	1 pulley	these are in a large box labelled "DYNAMIC CARTS AF22-2-13"
AF34-14	ruler	used to level the string
either AF35 (shelf) or AE82 (drawer)	1 metal ball (any size)	used to check the level of the track
AF44	larger weights	These are to ride the cart. Check with instructor: EITHER an assortment of 100-500 gram, cylindrical masses OR 2 black rectangular masses that fit into the cart (like the 1-wooden block above).
AF44 or AL32-5	tiny weights	These are to transfer between the cart and the basket. Check with instructor: EITHER an assortment of 7-10 very small masses (2-5 grams) OR 10 pennies.
-	Pasco	Computer
At the front for students to share		
S224	functioning digital scales	(please verify that these function and are set to metric)
AA41, AA42	Gravity Protractor	Used to level the track. (There are enough to give these to everybody, but it is not clear if there are enough pendulua strings and bobs for everybody.)
AF4-4	Pendulum bob to hang with the gravity-protractor	Used to level the track. (There are enough protractors to give these to everybody, but it is not clear if there are enough pendulua strings and bobs for everybody.)

Table 1.4.1: Equipment Needed: Newton's 2nd Law on a Linear Track with the Sonic Ranger




Figure 1.4.2: The equipment setup.

(Updated: October 12, 2017)

A digital version of the lab should be found at <http://physics.thomasmore.edu/PHY121Lab/c-Newton.html>

A PDF version might be found at [Newton.pdf](#)

1.5 Dry Sliding Friction

Most of this equipment is the same as [Newton's 2nd Law on a Linear Track with the Sonic Ranger](#), but it is not exactly the same.

Location	Equipment	Notes
For Each Lab Station		
S224	1 motion sensor (AKA “sonic ranger”)	Should have a black-yellow plug
S224	1 force transducer	Should be the thing shaped to sit on top of the wooden cart-block and it has a loop-screw at one side to tie a string to
AF12 (tube)	1 track	.
AF22-2-13	1 cart	these are in a large box labelled “DYNAMIC CARTS AF22-2-13”
AF22-2-13	1 wooden cart-block	these are in a large box labelled “DYNAMIC CARTS AF22-2-13”
AF22-2-13	light plastic bucket	These might already be attached to the string in a large box labelled “DYNAMIC CARTS AF22-2-13”
(attached to plastic bucket?)	string	There should be pre-cut string that is long enough to reach from the cart, over the pulley and to a hanging mass. About one meter long. This is likely in a large box labelled “DYNAMIC CARTS AF22-2-13”
AF22-2-13	1 pulley	these are in a large box labelled “DYNAMIC CARTS AF22-2-13”
AF34-14	ruler	used to level the string
either AF35 (shelf) or AE82 (drawer)	1 metal ball (any size)	used to check the level of the track
AF15	2-3 small wooden squares	used to level the track
AF44	larger weights	These are to ride the cart. Check with instructor: EITHER cylindrical masses (2 500-g, 4 200-g, 2 100-g) OR 2 black rectangular masses that fit into the cart (like the 1-wooden block above).
AF44 or AL32-5	tiny weights	These are to transfer between the cart and the basket. Check with instructor: EITHER an assortment of 7-10 very small masses (2-5 grams) OR 10 pennies.
AE21	string	This string is for pulling the force transducer, not the same as the string that is attached to the plastic bucket
AG15	1 thermal glove	Check with the instructor. This is only used as a soft glove to catch the cart in case it accelerates too quickly. Not all instructors want this.
-	Pasco	Computer
At the front for students to share		
S224	functioning digital scales	(please verify that these function and are set to metric)
AA41, AA42	Gravity Protractor	Used to level the track. (There are enough to give these to everybody, but it is not clear if there are enough pendula strings and bobs for everybody.)
AF4-4	Pendulum bob to hang with the gravity-protractor	Used to level the track. (There are enough protractors to give these to everybody, but it is not clear if there are enough pendula strings and bobs for everybody.)

Table 1.5.1: Equipment Needed: Dry Sliding Friction



(Updated: October 12, 2017)

A digital version of the lab should be found at <http://physics.thomasmore.edu/PHY121Lab/c-friction.html>

A PDF version might be found at [friction.pdf](#)

1.6 Centripetal Force

This describes the old equipment. New equipment (that plugs into the PASCO interface) has been purchased, but that equipment list has not been itemized yet.

Location	Equipment	Notes
For Each Lab Station		
AA13	1 caliper	<i>Ask instructor</i> if they want digital or analog (they might not care)
AA12	1 stop watch	Verify that they work and have good batteries. <i>You might need to purchase new batteries for department!!</i>
AD71	1 C-clamp	Important: There are small pieces of cardboard available to place between the clamp and the lab table so that you do not add dents to the table!
end of AD	1 thick metal rod	Mount horizontally off the table to hang the (spinning) carriage
AF14-1-6	1 centripetal force apparatus	(with card?)
AF46	1 0.5-kg hanger	.
AF46	Large cylindrical masses: 1 2-kg, 2 1-kg, 1 0.5-kg	.
AF44	Additional small cylindrical masses of various sizes	include at least one of each: 200-g, 100-g, 50-g, 20-g

Table 1.6.1: Equipment Needed: [Centripetal Force](#)



Figure 1.6.2: The equipment setup.

(Updated: October 12, 2017)

A digital version of the lab should be found at <http://physics.thomasmore.edu/PHY121Lab/c-centripetal.html>

A PDF version might be found at [centripetal.pdf](#)

1.7 Hooke's Law / Springs / Peer Review

Warning 1.7.1. This lab might be set-up by itself or it might be paired with [Pendulum / Peer Review](#). **Check with your instructor** to see if this is part of the “Peer Review” setup. If so, then discuss the setup. The instructor might want both Springs and Pendulum in both the morning and afternoon labs (in which case you need to verify how many setups of each lab and which tables), or the instructor might want one lab setup in the morning and a different lab setup in the afternoon.

Location	Equipment	Notes
For Each Lab Station		
AD41	1 2-slot-perpendicular clamp	This is used to connect the vertical thick rod to the horizontal small rod
AD71	1 C-clamp	Important: There are small pieces of cardboard available to place between the clamp and the lab table so that you do not add dents to the table!
AD82	1 clamp	This is used to hold the track in place. Note the size of rod that fits through the hole on this clamp. CAUTION: Do not over-tighten this clamp. <i>The track should not be bent by the clamp!</i>
end of AD	1 thick metal rod	Mount vertically off the table
end of AD	1 small rod	This is used to connect the clamp to the thick rod
AF12	1 track	This is used to measure the elongation of the spring.
AF44	hangers and various masses	each table should have a similar distribution of large and small mass sizes (Check with instructor: Your instructor might prefer to use the AF34 box of black hanging cylinders.)
At the front for students to share		
AF43	2 boxes of colored springs	Students will be allowed to select their spring. Check with instructor: Your instructor might want specific springs distributed to the tables

Table 1.7.2: Equipment Needed: [Hooke's Law / Springs / Peer Review](#)



(Updated: October 12, 2017)

A digital version of the lab should be found at <http://physics.thomasmore.edu/PHY121Lab/c-springs.html>


A PDF version might be found at [springs.pdf](#)

1.8 Pendulum / Peer Review

Warning 1.8.1. This lab might be set-up by itself or it might be paired with [Hooke's Law / Springs / Peer Review](#). **Check with your instructor** to see if this is part of the “Peer Review” setup. If so, then discuss the setup. The instructor might want both Springs and Pendulum in both the morning and afternoon labs (in which case you need to verify how many setups of each lab and which tables), or the instructor might want one lab setup in the morning and a different lab setup in the afternoon.

Location	Equipment	Notes
For Each Lab Station		
AA??, maybe AA41?	1 plastic protractor	To measure the amplitude of the pendulum swing.
AA12	1 stop watch	Verify that they work and have good batteries. <i>You might need to purchase new batteries for department!!</i>
AD21	1 3-notch black “arm” support	See picture below for how to mount the string inside the clamp. Do not just wrap it around the screw; clamp it inside the metal piece, which you should note has a nub holding it in place. (Take the metal square off, you'll see what I mean.)
AD71	1 C-clamp	Important: There are small pieces of cardboard available to place between the clamp and the lab table so that you do not add dents to the table!
end of AD	1 thick metal rod	Mount vertically off the table
AF44	hangers and various masses	each table should have a similar distribution of large and small mass sizes (Check with instructor: Your instructor might prefer to use the AF34 box of black hanging cylinders.)

Table 1.8.2: Equipment Needed: [Pendulum / Peer Review](#)



(Updated: October 12, 2017)

A digital version of the lab should be found at <http://physics.thomasmore.edu/PHY121Lab/c-pendulum.html>

A PDF version might be found at [pendulum.pdf](#)

1.9 Next Lab

Location	Equipment	Notes
For Each Lab Station		
AF12	1 track	.
At the front for students to share		
.		

Table 1.9.1: Equipment Needed: [Next Lab](#)

Figure 1.9.2: The equipment setup.

(Updated: October 12, 2017)

A digital version of the lab should be found at <http://physics.thomasmore.edu/PHY121Lab/c-labname.html>

A PDF version might be found at [labname.pdf](#)

1.10 Human Forearm

Location	Equipment	Notes
For Each Lab Station		
AD??	2 rod clamps (size depends on the rods you picked)	These connect the vertical rods to the horizontal rod
AE??	clamps for the meter-sticks	These are the little black things that slide onto a meter-stick and have a metal piece that a mass-hanger can hang off of.
End of AD	3 rods (2 long, 1 short)	Note that one rod goes in the table, so you need to notice if each table has a hole or a screw for this rod. The short rod is for stability.
End of AF	1 2-meter stick	Option You might use the rod with a base (and lead bricks at AG16) or a rod with a C-clamp (AD71), but the clamp might get in the way.
AF42	1 spring	(any color)
AF44	hangers and various masses	.
AF34	1 meter stick	(might need a 1/2-meter stick)

Table 1.10.1: Equipment Needed: [Human Forearm](#)




Figure 1.10.2: The equipment setup.

(Updated: February 27, 2013)

A digital version of the lab should be found at <http://physics.thomasmore.edu/PHY121Lab/c-equilibrium.html>

A PDF version might be found at [equilibrium.pdf](#)

Class 2

**PHY122L: General Physics
(algebra-based, spring)**